

## Additional information

### Patient safety

The instruments comply with IEC 60601-1:1988, IEC 60601-1-2:2001, IEC 60601-2-23:1999.

The following test house has approved the instrument: CSA in Canada according to CAN/CSA-C22.2 No. 601.1-M90, 601.151-94, 601.1B-98, and UL std. No. 601.2.23-98 and 60601-1.

Type BF equipment (body floating)



### CE

This product complies with the requirements of the Medical Device Directive 93/42/EEC June 1993 (Class IIa).

### EMC

This product complies with the requirements of the harmonized standard EN60601-1-2:2001, "Collateral standard: Electromagnetic compatibility – requirements and test".

### EMC emission

EN55011:1998, level A.

### EMC immunity

As stated in EN60601-1-2:2001, the immunity has been tested according to the IEC 61000 series.  
(See also IEC 60601-2-23:1999).

### Performance

This product complies with the IEC 60601-2-23:1999.

### Materials

All materials are latex-free.

### Languages

Danish, Dutch, English, French, German, Greek, Italian, Japanese, Portuguese, Russian, Spanish and Swedish.

## Make the most of transcutaneous monitoring

For more information on Radiometer's transcutaneous monitors and support services, visit [www.radiometer.com/tc](http://www.radiometer.com/tc).

For more clinical information on transcutaneous monitoring, visit [www.bloodgas.org](http://www.bloodgas.org).

## Sales companies:

<i>Country:</i>	<i>Radiometer representative:</i>	<i>Telephone:</i>	<i>E-mail:</i>	<i>Website:</i>
Australia:	Radiometer Pacific Pty. Ltd.	(+61) 3 9259 2222	<a href="mailto:sales@radiometer.com.au">sales@radiometer.com.au</a>	<a href="http://www.radiometer.com.au">www.radiometer.com.au</a>
Canada:	Radiometer Canada	(+1) 877 414 0447 (toll free)	<a href="mailto:info@radiometercanada.com">info@radiometercanada.com</a>	<a href="http://www.radiometer.ca">www.radiometer.ca</a>
China:	Radiometer China	(+86) 21 6128 6466	<a href="mailto:pol@radiometer.com.cn">pol@radiometer.com.cn</a>	<a href="http://www.radiometer.cn">www.radiometer.cn</a>
Denmark:	Radiometer Danmark	(+45) 38 27 28 29	<a href="mailto:rdan@radiometer.dk">rdan@radiometer.dk</a>	<a href="http://www.radiometer.dk">www.radiometer.dk</a>
France:	Radiometer S.A.S.	(+33) 1 49 44 35 50	<a href="mailto:marketing@radiometer.fr">marketing@radiometer.fr</a>	<a href="http://www.radiometer.fr">www.radiometer.fr</a>
Germany:	Radiometer GmbH	(+49) 2154 8180	<a href="mailto:info@radiometer.de">info@radiometer.de</a>	<a href="http://www.radiometer.de">www.radiometer.de</a>
Ireland:	Radiometer Ireland Ltd.	(+353) 1 888 3611	<a href="mailto:sales@radiometer.ie">sales@radiometer.ie</a>	<a href="http://www.radiometer.ie">www.radiometer.ie</a>
Japan:	Radiometer K.K.	(+81) 3 5777 3500	<a href="mailto:salesdep@radiometer.co.jp">salesdep@radiometer.co.jp</a>	<a href="http://www.radiometer.co.jp">www.radiometer.co.jp</a>
The Netherlands:	Radiometer Nederland BV	(+31) 79 361 4593	<a href="mailto:info@radiometer.nl">info@radiometer.nl</a>	<a href="http://www.radiometer.nl">www.radiometer.nl</a>
New Zealand:	Radiometer New Zealand	(+64) 9 574 1400	<a href="mailto:sales@radiometer.co.nz">sales@radiometer.co.nz</a>	<a href="http://www.radiometer.co.nz">www.radiometer.co.nz</a>
Poland:	Radiometer Sp. z o. o.	(+48) 22 518 02 40	<a href="mailto:info@radiometer.pl">info@radiometer.pl</a>	<a href="http://www.radiometer.pl">www.radiometer.pl</a>
Portugal:	Radiometer Ibérica, S.A.	(+351) 214 12 39 70	<a href="mailto:rapor@radiometer.es">rapor@radiometer.es</a>	<a href="http://www.radiometer.pt">www.radiometer.pt</a>
Spain:	Radiometer Ibérica, S.A.	(+34) 91 655 99 50	<a href="mailto:resp@radiometer.es">resp@radiometer.es</a>	<a href="http://www.radiometer.es">www.radiometer.es</a>
Switzerland:	Radiometer GmbH	(+41) 44 723 38 60	<a href="mailto:info@radiometer.ch">info@radiometer.ch</a>	<a href="http://www.radiometer.ch">www.radiometer.ch</a>
United Kingdom:	Radiometer Ltd.	(+44) 1293 517 599	<a href="mailto:sales@radiometer.co.uk">sales@radiometer.co.uk</a>	<a href="http://www.radiometer.com">www.radiometer.com</a>
USA:	Radiometer America Inc.	(+1) 800 736 0600 (toll free)	<a href="mailto:info@radiometeramerica.com">info@radiometeramerica.com</a>	<a href="http://www.radiometeramerica.com">www.radiometeramerica.com</a>
Other countries:	Radiometer International Sales Division	(+45) 38 27 38 27	<a href="mailto:isd@radiometer.dk">isd@radiometer.dk</a>	<a href="http://www.radiometer.com">www.radiometer.com</a>

Data subject to change without notice.

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, PICO and CLINITUBES are trademarks of Radiometer Medical ApS.

# TCM400 Specifications



# Parameter configuration

Type	Parameters	Units	Ranges
Measured			
Transcutaneous oxygen tension	$tcpO_2$	mmHg kPa	0-2000 0-266.7
Electrode heating	Power	mW	10-650
Calculated			
Regional Perfusion Index	RPI	-	0-3

## Monitor data

Display options:	Normal view (numeric), Trend table view, Trend curve view
Display update:	Every 2 sec.
Print reports:	Trend table, Trend curve
Barometer:	375-825 mmHg, 50-110 kPa
Calibration:	Ambient air or calibration gas
Start-up time:	Max. 1 min.
Time:	Date, Clock
Timer:	Range: 0-99 hours, increments: 1 sec/1 min

### Power requirements

Voltage:	90-264 VAC
Frequency range:	47-63 Hz
Power consumption:	70 VA (max.)

### Monitor battery

Type:	Rechargeable Pb battery
Duty period:	1 hour typical per charge at 25 °C
Recharging time:	Approx. 8 hours at 25 °C

### Sensor temperature settings

Setting:	37-45 °C
Increments:	0.5 °C
Accuracy:	± 0.1 °C (excluding sensor)



## Dimensions

Monitor and module dimensions					Electrode dimensions (E5250)		
	Monitor		O <sub>2</sub> Module				
Height:	16 cm	6.3"	3.5 cm	1.4"	Diameter:	Ø 15 mm	0.6"
Width:	30.8 cm	12.1"	14.5 cm	5.7"	Height:	11.3 mm	0.44"
Depth:	23 cm	9.1"	14.8 cm	5.8"	Weight:	2.9 g	0.1 oz
Weight:	4 kg	8.8 lbs	0.22 kg	0.5 lbs			

## IT solution

### Computer specifications

6½" TFT, full VGA (640 x 480) color touch screen  
 AMD Geode GX1, 300 MHz (Pentium Class)  
 Windows CE 4.2  
 64 MB RAM  
 48 hours' storage of measuring data

### Interface possibilities

Serial line: EIA232, (RS232)  
 Printer output: USB 1.0, HP PCL3 or higher

## Sensor specifications

### E5250

O<sub>2</sub> cathode: 25 µm Platinum  
 O<sub>2</sub> anode (reference): Silver  
 Measuring principle: Transcutaneous Clark-type O<sub>2</sub> electrode

## Sensor performance (at 43 °C)

### E5250

Response time pO<sub>2</sub>: < 11 sec (10 % to 90 % response)  
 Max. drift per hour: < 1.0 %  
 Linearity at 0 % O<sub>2</sub>: Better than 1 mmHg  
 Linearity at 90 % O<sub>2</sub>: Better than 25 mmHg (equal to 4 %)

## Accessories

### Fixation rings (904-891)

Diameter: 30 mm  
 Adhesive material: Medical grade Acrylic adhesive  
 Ring material: PVC  
 Contact solution: 1.2-propanediol and deionized water

### Membranes (904-308)

Membrane material: PP membranes  
 Electrolyte solution: 1.2-propanediol, potassium chloride, sodium hydrogen carbonate and deionized water

### TCC3 calibration unit (optional)

Height: 8 cm 3.15"  
 Width: 12 cm 4.7"  
 Depth: 23 cm 9.0"  
 Weight: 1.9 kg 4.2 lbs

Calibration: 1-point  
 Calibration gas: 5 % CO<sub>2</sub> and 20.9 % O<sub>2</sub>, balance N<sub>2</sub>

Gas flow: 8 mL/min ± 2 mL/min  
 Automatic shut-off: After 5, 10, 15, 20, 50 minutes, as desired

Battery type: 1.5 Alkaline, IEC type LR6

