



Vitalograph COMPACT

Sophisticated, integrated cardio-respiratory

diagnostic workstation

The Vitalograph COMPACT™ Expert is a medical workstation that integrates high quality spirometry, 12-lead ECG, Pulse Oximetry, COPD assessment, blood pressure measurement, challenge testing and more in a flexible, integrated solution. COMPACT features secure data management using data encryption in

Vitalograph is a world leading provider of outstanding quality cardio-respiratory diagnostic devices, clinical trial services and medical equipment servicing. With a pioneering heritage of excellence spanning half a century Vitalograph continues to make valuable contributions to effective medical care and enhanced quality of life.


COMPACT Workstation

SOL databases.

The COMPACT Expert workstation integrates the management of your subjects, visits and all of your clinical testing. Functionality can be customised for clinical studies or other major projects. The Vitalograph COMPACT offers:

- Networking of multiple COMPACT devices in their own network or to a central database on your existing networks
- Simple icon driven menu on a large colour touch screen with on-screen keypad and optional biometric ID for users and patients. A standard mouse and keyboard may be used if preferred
- Intuitive user interface with workflow orientated presentation in four primary information panels: subjects, test sessions, information and session test panel
- Automatic storage and traceability of all test data and events, including calibration logs and audit trails
- Designed for secure encrypted data transmission via the internet to Vitalograph hosted central servers for reporting, archiving, analysis, review, over-reading, etc.

- Manual data entry and data export in PDF or CSV format
- Automated, secure data management designed to comply with the technical requirements of FDA 21 CFR 11, 95/46/EC and HIPAA compliance










Spirometry

The COMPACT Expert is the latest evolution from Vitalograph's range of world leading spirometers. With vastly increased power, and storage capacity to accommodate many other types of medical device connectivity, the integrity of the spirometry device is not lost but enhanced. The spirometry module features:


- Scalable volume/time & flow/volume curves using single-breath or multi-breath spirometry
- Automatic and scalable FEV₁ trend chart view as subject is selected prior to testing, designed to reflect the trend of FEV₁ over the lifetime of the subject
- Optional choice of child incentives with sound effects and transparency so that real-time curves are simultaneously viewable
- Pre-post spirometry comparison for bronchodilator or longer-acting treatments
- User selection from a choice of over 50 test parameters

Vitalograph COMPACT for Communication and **Training**

- Ethernet x 2, Serial ports x 2
- USB 2.0 x 4
- Bluetooth® and WIFI Dongle option
- Audio out (stereo) and microphone input
- DVI video output for connection to a monitor or projector



User preference for comparison with a 'normal' population in % predicted or Z-score (SDS)



Vitalograph COMPACT High Quality Flexible Reports

On any network or USB printer







The COMPACT Expert with ECG offers full 12-lead electrocardiography with very low occurrence of lead artefact. Analysis and interpretation come as standard using the world class University of Glasgow interpretation algorithm which provides detailed diagnostic suggestions and evaluation of results. ECG measurements are captured and transmitted wirelessly for display in real-time and data is automatically saved for recall, reporting or printing. Transmitting ECG traces wirelessly is fast and safe with no leads to trip over or create artefacts from movement or electrical sources. The built-in display provides information on electrode contact problems, heart rate, battery status, communication status and digital pacemaker detection. The COMPACT Expert with ECG features:



- Operator controlled selection of a 10-s sample trace, but can capture and store several minutes of 12-lead ECG recording
- Visual electrode contact indication to assist operator
- Wireless data transmission via Bluetooth (range: approx. 10m)
- High quality multi-channel and 12-lead overview for viewing and printing
- Print configurable 1 or 2 page reports
- View traces as 1 x 12; 2 x 6 or 3 x 4 plus rhythm strip
- Calculation of the heart rate with acoustic and visual signal
- Full range of parameters are available (e.g. P duration, PQ interval, QRS period)

- Automatic pacemaker detection
- Rhythm analysis







Bronchial Challenge (ATS and Mannitol)

The causes of asthma and allergy can be difficult to diagnose, but bronchial provocation with specific sensitisers, methacholine or mannitol is of value to identify the extent and cause of bronchial hyper-responsiveness.



Pulse Oximetry


The Spirotrac® Pulse Oximetry Software feature allows you to plug an SpO2 sensor directly in to the USB port on your COMPACT Expert to deliver accurate and high quality oximetry alongside spirometry.

This calculates, trends and reports all primary parameters including minimum, maximum and median SpO2/pulse rate. Capable of recording quick assessments or longer term recordings and up to four Tx values can be selected by the user.



COPD Assessment

- Spirometry and pulse oximetry data is automatically or manually integrated using the ATS/ERS 6-minute walk test protocol
- 6-minute walk test Wizard runs in 10 logical stages
- The BODE* Index (*Body Mass Index, Obstructive Index, Dyspnoea Score, Exercise Tolerance) predicts COPD mortality based on FEV₁, 6MWT, BMI, MMRC and BORG Scale





Blood Pressure

Subject BP data can be transferred to the database from a Bluetooth Blood Pressure device or entered manually.



Weighing Scale

Subject weight data can be transferred to the database from a Bluetooth Medical Weighing Scales or entered manually.

italograph

UK Sales Vitalograph Ltd. **Maids Moreton Buckingham** MK18 1SW England Tel: 01280 827110 Fax: 01280 823302 e-mail: sales@vitalograph.co.uk www.vitalograph.co.uk

Export Sales Vitalograph Ltd. Maids Moreton Buckingham MK18 1SW England Tel: +44 1280 827120 Fax: +44 1280 823302 e-mail: sales@vitalograph.co.uk e-mail: info@vitalograph.de www.vitalograph.eu

Vitalograph GmbH Rellinger Straße 64a 20257 Hamburg Germany Tel: (040) 547391-0 Fax: (040) 547391-40 www.vitalograph.de

Vitalograph Inc. 13310 West 99th Street Lenexa Kansas 66215 USA Toll Free: 800 255 6626 Tel: (913) 730 3200 Fax: (913) 730 3232 e-mail: vitcs@vitalograph.com www.vitalograph.com

Vitalograph (Ireland) Ltd. **Gort Road Business Park Fnnis** Co Clare Ireland Tel: (065) 6864100 Fax: (065) 6829289 e-mail: sales@vitalograph.ie www.vitalograph.ie

www.vitalograph.co.uk



Technical Specification:

Product: COMPACT Model No: 6600

Net Size: 375mm x 235mm x 110mm **Gross Size:** 461mm x 297mm x 226mm Weight: Net: 2.5kg; Gross: 4.3kg

Spirometry:

Flow Detection Principle: Fleisch type

pneumotachograph

Spirometry Parameters: Selectable, in excess of

Spirometry Predicted Values: Selectable

Spirometry Test Types: Single breath tests, flow/ volume loops, multi-breath testing, tidal breathing

Volume Accuracy: Better than +/ 3%

Recommended Operating Temperature Range:

17°C - 37°C

Spirometer Performance Standards:

ISO26782:2009; ISO 23747:2009;

ATS/ERS 2005

ECG:

BT12: 12-channel, 3 leads acc. to Einthoven, 3 acc.

to Goldberger, 6 Wilson leads

Power supply: 2 x AA batteries (1.5 V) or

rechargeable batteries (1.2 V)

Runtime when used constantly approx. 14 h

Electrodes: Standard ECG electrodes

Operating Temperature Range:

10°C - 40°C

ECG Safety & Performance Standard:

IEC 60601-2-25

Data Storage: Stores up to 10,000 subjects

Input Voltage/Frequency: 110-250VAC 50/60Hz

PowerSAFE Output Voltage: 19VDC **Storage Humidity:** 10% to 95% **Storage Temperature:** 0°C to 50°C Safety Standards: IEC 60601-1:2005 **Medical Safety Standard:** Medical Devices Directive 93/42/EEC1993 (as amended)

Class: Ila MDD

Designed & Manufactured Under:

ISO 13485:2003, FDA 21CFR820 & CMDR

Communications: USB 2.0 x 4, Serial RS232 x 2,

Ethernet x2

DVI video display interface

Ordering Info:

67800 6600 Vitalograph COMPACT with Fleisch Pneumotachograph includes Software modules for Spirometry, ECG, SpO2, BP, Weight, Bronchial Challenge Testing,

67900 6600 COMPACT with ECG

COMPACT Compatible Devices:

70400 Pulse Oximetry Finger Probe USB

70401 Biometric Scanner USB

70403 Blood Pressure Monitor Bluetooth

70404 Medical Weighing Scales Bluetooth

Other Vitalograph medical devices that communicate with COMPACT:

79400 2120 Hand Held In2itive™ Spirometer

65500 6000 ALPHA Touch™

40400 4000 Respiratory Monitor asma-1[™] USB

Consumables and Accessories

41301 Disposable ECG Electrodes (30)

28350 2820 BVF™ Bacterial/Viral Filters (50)

20242 2024 SafeTway® mouthpieces (200)

20980 2024 Mini SafeTway mouthpieces (50)

20303 2030 Noseclip disposable (200)

36020 2040 Precision Syringe 3-L 41302 ECG Simulator





